

Puryx™ Rapid

Simply rapid DNA / RNA sample purification

Puryx™ Rapid Extraction Kits

The straight-forward protocol of Puryx™ Rapid has a low number of steps to enable sample purification in a fraction of the time.


All reagents required to perform the extraction are included in the box.

Perform manually or with open semi-automated robotic platforms.



Puryx™ Rapid sample purification process

- Purify within just 25 minutes
- Just one step for lysis and binding and a single wash step
- No need for additional reagent preparation
- DNA and RNA extraction
- Nasal, pharyngeal, and throat swabs



Ordering

93620	Puryx™ Rapid DNA / RNA Extraction Kit
93621	Deep well plate 1.1 mL (50)
93611	Deep well plate 2.2 mL (100)
Enquire	Universal magnetic plate
Enquire	300 mL reservoirs, sterile

Contact us now to find out how Puryx™ Rapid will help you to achieve faster sample purification.

©2020, AusDiagnostics Pty Ltd. Puryx is a trademark of AusDiagnostics Pty Ltd. Puryx™ Rapid extraction kits are intended for use by trained professionals to purify nucleic acids from clinical samples, and are for *in vitro* diagnostics use. FPURYXRPD_01; 2020 OCT 12.



www.ausdiagnostics.com

AusDiagnostics Pty Ltd
290–292 Coward St
Mascot NSW 2020
Australia
T: +61 2 9698 8030

AusDiagnostics UK Ltd
Unit 3, Anglo Business Park
Asheridge Rd, Chesham,
Bucks, HP5 2QA, UK
T: +44 0 1494 300 121

AusDiagnostics NZ Ltd
7C Douglas Alexander Pde
Rosedale AK 0632
New Zealand
T: +64 9 886 9717

AusDiagnostics USA Inc
Bldg 2, 1367 Brumlow Ave
Southlake TX 76092
USA
T: +1 940 279 2171