

20th June 2017

Advisory Notice: RSV background update

There is an update to a product you order. AusDiagnostics routine monitoring of results files resulted in an update to the template for Respiratory 8-plex panels. This update will reduce the excessive background RSV 'check' calls.

Please install the new template of the product you use from the AusDiagnostics website (www.ausdiagnostics.com) or ask your sales representative.

Please find the products affected below:

Name of IVD products affected:

- REF 20081 VER 05 - Influenza + RSV (8-well)
- REF 20614 VER 02 - Respiratory Pathogens (8-well)

Please don't hesitate to contact us on 02 96988030 or at support@ausdiagnostics.com if you have any further queries or require any assistance with the update procedure.



AusDiagnostics Pty Ltd
205 Victoria Street
BEACONSFIELD
NSW 2015 Australia
Phone: +612 9698 8030
Fax: +612 9319 1661



AusDiagnostics UK Ltd
Unit 3
Anglo Business Park, Asheridge Rd,
Chesham, Bucks, HP5 2QA UK
Tel: +44 (0) 1494 300121
Mob: +44 (0) 7896 534843