

20 December 2016

## Advisory Notice – Faecal Pathogens B (16-well)

There will be a future update to a product you order, a new IFU will be released to reflect these changes. Please find the product, and details of the change below:

### Faecal Pathogens B (16-well)

**Previous catalogue number of product:** REF 25033 VER 03

**New catalogue number of product:** REF 25033 VER 04

**New Document ID of IFU:** 25033r02

**Regulatory status of updated product:** IVD use in Australia and Europe

#### Reason/s for change:

- The Salmonella has been updated to reduce the incidence of false positives. The gene target is unchanged.
- Noro 1 and Noro 2 assays updated to remove background. The gene target is unchanged.
- The RT conditions have been updated to those used in Respiratory and CSF panels. The gene target is unchanged.

#### Date effective of change: 3rd January 2017

Please don't hesitate to contact us in Australia on +61 (0) 2 96988030 (Australia, New Zealand, Oceania & Asia) or in the UK on +44 (0) 1494 300121 if you have any further queries or require any assistance with the updated products.