

21 December 2016

Advisory Notice – Format of *Easy-Plex™* 96 products

This note is to let you know that AusDiagnostics will adopt the following packaging and prices for *Easy-Plex™* 96 products effective 01 January 2017. This is an interim arrangement while specialised plastic ware is manufactured.

16 well kits:

Step 1 tubes will be sold in bags of six 8-well strips sufficient for 48 samples. These will be sealed with caps and placed in a resealable bag. This is sufficient for 6 months use at two samples per week. Only remove as many step 1 tubes as you need (the 8-well strips may be cut into individual tubes). Do not uncap tubes that you are not using. Shelf life is 6 months if handled in this way. 48 dilution tubes will be supplied with each step 1 box. These are used with the 24 x 2 ml loading block in the *Easy-Plex™* processor.

Step 2 boxes will contain eight 96-well plates sufficient for 48 samples. Each plate can be used for 6 samples but must be used once the pack is opened even if only one sample is run. The price of the plates has been reduced so that the price per sample is similar for 2 or more samples.

This kit requires three demi master mixes to run 6 samples

24 well kits:

Step 1 tubes will be sold in bags of six 8-well strips sufficient for 48 samples as for 16-well kits.

Step 2 boxes will contain twelve 96-well plates sufficient for 48 samples. Each plate can be used for 4 samples but must be used once the pack is opened even if only one sample is run. The price of the plates has been reduced so that the price per sample is similar for 2 or more samples.

This kit requires two demi master mixes to run 4 samples

The new IFU 9350r04 will be available early January 2017

The new versions of kits with new templates will be released to support the described changes.

The price per sample for *Easy-Plex™* 96 is compared to High-Plex in this table.

Samples	16-plex	EP 16-well	EP 24-well
Box \$	\$4950.00	\$1416.00	\$1420.00
1	\$80.83	\$47.00	\$52.00
2	\$48.33	\$33.50	\$38.50
3	\$37.50	\$34.00	\$39.00
4	\$32.08	\$30.50	\$35.50
5	\$28.83	\$31.40	
6	\$26.67	\$29.50	

Please don't hesitate to contact us at support@ausdiagnostics.com if you require further information. Thank you for your continuing support.