

1st March 2017

Advisory Notice: Melt range update – HSV-1

There is an update to a product you order, due to AusDiagnostics routine monitoring of melt temperature resulted in an update to the melt temperature range of the HSV-1 assay.

Please install the new template of the product you use from the AusDiagnostics website (www.ausdiagnostics.com) or ask your sales representative.

Please find the products affected below:

Name of IVD products affected:

- REF 27091 VER 08 – Herpes, Enterovirus and Adenovirus (8-well)
- REF 27093 VER 01 – Viral (8-well)
- REF 27092 VER 01 – Viral CSF (8-well)

- REF 67091 VER 02 – Herpes (6-well)
- REF 67060 VER 09 – CSF (14-well)

Name of RUO products affected:

2705002, 2711201, 3760104, 3714101, 3705002, 3719003, 3708002, 3702503, 3702105, 3702305, 7705001, 7702401, 8709101

Please don't hesitate to contact us on 02 96988030 or at support@ausdiagnostics.com if you have any further queries or require any assistance with the update procedure.



AusDiagnostics Pty Ltd
205 Victoria Street
BEACONSFIELD
NSW 2015 Australia
Phone: +612 9698 8030
Fax: +612 9319 1661



AusDiagnostics UK Ltd
Unit 3
Anglo Business Park, Asheridge Rd,
Chesham, Bucks, HP5 2QA UK
Tel: +44 (0) 1494 300121
Mob: +44 (0) 7896 534843