

27th of July 2017

Advisory Notice: Update to Synthetic Positive Control for Faecal Products

Dear Valued Customer,

We will soon be releasing changes to our range of panels for the detection of Faecal Pathogens. These changes include an update to our campylobacter assay so that it is specific for pathogenic *Campylobacter* species only.

As a result of this, we will be updating the Synthetic Positive Control for Faecal Panels.

If you currently use one of the following panels, please be advised that you should only use Synthetic Positive Control for Faecal Panels (REF 91031) with an expiry date of February 2018 onwards.

<p>Faecal Pathogens A (16-well) Catalogue number of product: REF 25031 (VER 06)</p>	<p>Faecal Pathogens M (16-well) Catalogue number of product: REF 25039 (VER 06)</p>
<p>Faecal Pathogens B (16-well) Catalogue number of product: REF 25033 (VER 05)</p>	<p>Faecal Bacteria and Parasites (8-well) Catalogue number of product: REF 25038 (VER 03)</p>
<p>Faecal Pathogens C (16-well) Catalogue number of product: REF 25036 (VER 06)</p>	

If you currently have stock of this product with an expiry date that falls between August 2017 and February 2018, please contact AusDiagnostics for a free of charge replacement. If you are unsure as to whether you have the correct product, please do not hesitate to contact us at support@ausdiagnostics.com or through the details below.

Effective date of change: 28th of July 2017.

AusDiagnostics Pty Ltd

205 Victoria Street
BEACONSFIELD
NSW 2015 Australia
Phone: +612 9698 8030
Fax: +612 9319 1661

AusDiagnostics UK Ltd

Unit 3
Anglo Business Park, Asheridge Rd,
Chesham, Bucks, HP5 2QA UK
Tel: +44 (0) 1494 300121
Mob: +44 (0) 7896 534843

AusDiagnostics Pty Ltd

205 Victoria Street, Beaconsfield NSW 2015 Australia ABN 84 612 059 686
P +61 2 9698 8030 F +61 2 9319 1661 E support@ausdiagnostics.com ausdiagnostics.com