

Comprehensive Diagnosis for Enteric Pathogens

- ❑ Covers seven of the most common bacterial pathogens and Shiga toxins
- ❑ Three obligatory pathogenic parasites
- ❑ Five families of enteric viruses
- ❑ All pathogens detected from a single sample
- ❑ Ready-to-use reagents
- ❑ Automated process and results calling
- ❑ Barcode traceability and LIMS compatibility
- ❑ Minimum hands-on time

Tandem-Plex™ Faecal Pathogens M

Simultaneous detection of 18 pathogen targets from a single sample.



PANEL TARGETS	
Bacteria Salmonella spp. Shigella spp. Shigatoxin 1 & 2 Campylobacter jejuni & coli Yersinia enterocolitica & pseudotuberculosis Aeromonas spp. Clostridium difficile	Parasites Entamoeba histolytica Cryptosporidium hominis & parvum Giardia duodenalis (G.lambliia)
	Viruses Norovirus 1 / Norovirus 2 Rotavirus / Astrovirus Adenovirus / Sapovirus

CONTROLS
Artificial sequence for internal control of amplification
Extraction control*
Synthetic positive control*

* optional

E.coli O157 target available in a similar product REF.25033 Faecal Pathogens B